



PROTEIN, BLOOD & STAIN REMOVER

B-Gone Blood & Stain

Protein, Blood, and Stain Remover

CORE 3

Core 3: Personal Attention

Every customer requires special attention and has their own demands put on them. Our staff will design and implement a program designed to fit your needs.

Core 3: Experience

With our full staff of formulators, chemists, & marketing professionals we have been bringing cutting edge products to market for over 70 years.

Core 3: Technology

Through daily contact with industry professionals, we are continuously updating and improving our products to make sure you perform at your best.

Synthetic Labs offers industry-leading solutions for businesses and professionals who demand exceptional results.



PRE-SPOTTER

B-Gone Blood & Stain

Protein, Blood, and Stain Remover

What is B-Gone Blood & Stain?

A powerful laundry spotter that has been specially formulated with emulsifiers, alkali builders, and solvents to quickly penetrate, breakdown and remove blood & protein stains.

Benefits:

- Loosens tough embedded oils.
- Free rinsing formula.
- Extremely versatile product.

Precautionary Statements:

For institutional or commercial use.

Spotting Board Instructions:

1. Identify the stain
2. Blot the stain gently with a clean cloth to absorb excess material.
3. Apply stain remover directly to the stain.
4. Let the stain remover sit for the recommended time (usually 15-20 minutes). The longer, the better
5. Gently rub the fabric together or use a soft brush to work in the remover.

Wash in hot water, 140 for colors and 165 for whites..

Dry Cleaning Instructions:

Apply directly to stain.

Launder through a dry cleaning machine.

Tech Specs:

Color.....Clear
 Fragrance.....Ammonia
 Form.....Liquid
 Phosphates.....None
 Flash Point.....212°F
 pH.....9
 Specific Gravity...1.04 g/m³
 Viscosity.....Moderate
 Foam.....Moderate
 Biodegradable.....100%
 Shelf Life.....1 Year

Available Packages

9001-6QT

6x32 OZ. Quart Case



Toll Free: 800-255-4050

24 Victory Lane, Dracut, MA 01826

www.syntecpro.com